



Display  
  
12.1"

5  
Basic  
Parameters

Model  
ECO

User  
Friendly  


  
USB  
Connect

*“Comprehensive and Powerful  
Monitoring Solutions”*

# *Multi-Para Patient Monitor*

# Vital Signs Show Perfectly



**Trend Graph**



**Big Font Interface**



**OXI-CRG Display**



**Durable Probes  
Connector with High  
Quality**



**Max7 Wave form display  
simultaneously**

- 12.1" TFT Color Display
- Neonatal and Adult Usage
- 6 or more waveform along with numerical parameters
- Basic 5 Parameters ECG, SPO2, NIBP, RESP, TEMP.
- Facilitates adjustment & setting display Alarm.
- Display battery status indicator

## ECG

- Heart rate range from 15-120 bpm
- Auto, ECG or Plethysmograph as source for heart rate
- Sweep speed of 12.5/25/50mm/sec
- Display audio, visual alarm
- Display a range of 0-120 breath pm for respiration

## NIBP

- Based upon Oscillometric Method
- Modes: Manual/ Auto/ Stats
- Display systolic, diastolic and mean arterial pressure
- Systolic BP range from 0-300; Diastolic BP range from 10-200 and MAP from 20-220 mmHg

## TEMP

- Temperature range of 25-45 degree C
- Resolution of 0.1 degree C
- Accuracy of +/-0.2 degree C

## PULSEOXIMETRY

- Saturation range of 0-100% for adult child and neonate
- PR range of 30-250 bpm
- Plethysmographic waveform display

## TRENDS

- 48 hrs trend storage of all parameters in graphic and tabular
- Trend Graph, Big font interface, OXI-CRG display

## OPTIONAL

Dual IBP, ETCO2, Printer

## ACER BIOMEDICALS LLP

Office: Industrial Unit-80, Block-G, Sector-63, Noida, UP - 201301 (India)  
Ph.: +91-120-4166921, E-mail: info@acerbiomed.com  
Web.: www.acerbiomedicals.com, www.acerbiomed.com

Service Helpline  
18008894232

