

Multi-Para Patient Monitor




Model: CM-50

12.1" Inch



Professional	Accurate	Reliable	Smart
Simple	Cost-effective	Portable	Connectivity

Software interface diagram




Standard Interface



Large Font Interface



Trending Interface



Multi-guides with Screen Interface



Quad Split Screen Interface



NIBP

Method- Pulse wave oscillation method

Measure mode- Manual/automatic/continuous

Respiration Rate

Measurement Range

Adult	0 ~ 120brpm
Neonatal / Pediatric	0 ~ 150brpm
Resolution	1brpm
Accuracy	06brpm: unspecified 7~150brpm: ±2brpm

Heart rate

Measure range

Adult	15 ~ 300 bpm
Neonatal/Pediatric	15 ~ 350 bpm
Accuracy	1% or 1bpm, Take the head
resolution	1 bpm

Monitor the size and weight

Size	305×290×145 mm
weight	2.80Kg

Power Supply

External power requirements: a.c 100~240 V, 50/60 Hz, Maximum input power: 90VA , Fuse : T1.5AL 250V

Battery- 2.2 mAh 14.8V rechargeable battery

SpO₂

Measurement Range	0 ~ 100 %
Resolution	1 %
Accuracy	70% ~ 100% : 2 % <69% :unspecified

Pulse Rate

Measure and Alarm Range	20~250bpm
Resolution	1bpm
Accuracy	3bpm

Temperature

Sensor type	10 K series, 2.25 K series
Channel	2
Measure and Alarm Range	0 ~ 50 C
Resolution	0.1C
Accuracy	± 0.1°C(Do not include the sensor error)

Pulse Rate Range	40 ~ 240 bpm
-------------------------	--------------

ACER BIOMEDICALS LLP

Office: Industrial unit -80, Block-G, Sector-63, Noida, U.P.(India)- 201301
 +91-120-4166921, E-mail: info@acerbiomed.com
 Web: www.acerbiomedicals.com, www.acerbiomed.com

Toll-free No: 18008894232