



HbA1c Analyser HB-100


- Small and portable analyzer, easy to use.
- Tiny sample size 10 μ L, fingerstick sample and no pretreatment.
- Fast and accurate result in 2 minutes.
- Support blood glucose test, more perfect diabetes care.

CE

HbA1c Analyser
HB-100
 A Point of Care Analyser



3 SIMPLE STEPS



SPECIFICATIONS

Advantaged Boronate affinity chromatography for HbA1c test
Sample type : whole blood from finger tip or venous blood samples
Sample volume : 10pL
Measuring range : 4%-16% HbA1c
Imprecision : < 3%
Stability : 4-30°C (39.2-86 T); 18months validity period and 30 days on-board reagent stability
Report form : IFCC (mmol/mol) and NGSP/DCCT (%) values

Smart and easy-to-use analyzer

Dimensions (L*W*H) : 135mmx68mmx25mm
Weight : 0.2kg
Power supply : External power and integrated rechargeable lithium-ion battery (800mAh, standby time 48hours)
External Input/Output : RS232C (ASTM), Bluetooth
Memory : 160 results
Printer : Optional built-in thermal printer
Operating environment : 10-40 C (50-104 T), 30-70% R.H.(Without condensation)

Supporting Kit : Firstview HbA1c Quantitative Test (25 test kit)

Service Helpline
18008894232



Marketed By
Acer Biomedicals
 Office: Industrial Unit-80, Block-G, Sector-63,Noida UP- 201301, India ;
 Ph: +91-120-4166921,Email- info@acerbiomed.com; sales@acerbiomed.com;
 Web: www.acerbiomedicals.com ; www.acerbiomed.com