



PORTABLE SUCTION MACHINE 7E-H1



TECHNICAL SPECIFICATIONS		
Input power	:	90 VA
Max Vacuum	:	≥ 0.08 MPa
Max Air Flow	:	28 \pm 2 L/min (without jar) 22 \pm 2 L/min (with jar)
Operation noise	:	Approx. 60dB(A) (without vacuum)
Liquid Collection Jar	:	1000mL, 1piece
Operating time	:	30 minutes on, 30 minutes off

ACER BIOMEDICALS LLP

Office: Industrial unit -80, Block-G, Sector-63, Noida, U.P.(India)
+91 120 4166921, E-mail: info@acerbiomed.com
Web: www.acerbiomedicals.com, www.acerbiomed.com
Toll free: 18008894232

Service Helpline
+91 120 4166921

