



# Thermometer TH-C11



Non-contact  
measurement



Fever warning



Memorize data



Automatic shut-down



Quick temperature  
measurement

Product name	Infrared thermometer
Model No	TH-C11
Power supply	DC3V(2*AAA Cell)
Size	195*125*52MM
Resolution	0.1C/32.18F
Measuring range	34C(93.2°F)-43 C(109.4 F) (Body temperature mode)
Warranty	5years
Accuracy	+/- 0.2°C/32.36°F(within 35°C(95°F)~42°C(107.6°F)

## ACER BIOMEDICALS LLP

Office: Industrial unit -80, Block-G, Sector-63, Noida, U.P.(India)- 201301

+91-120-4166921, E-mail: info@acerbiomed.com

Web: www.acerbiomedicals.com, www.acerbiomed.com

Toll-free No: 18008894232