



Finger Tip Pulse Oximeter

FP-22




Accurate




Small & Elegance



Reliable & Stable



Weaving & Comfortable



Affinity & Design

Blood Oxygen Saturation Value

Pulse Rate



Blood Oxygen Volume Waveform

Histogram of Pulse Intensity

Main Specifications:

Power supply: 2AAA battery

Display: color OLED

Working current: 30mA

Resolution: +-2%

Dimensions: 59*36*34mm

Pulse rate(PR): 30-250 bpm

Blood oxygen saturation: 35-100%

Main Features:

- Easy Operation.
- Convenient detection with universal display.
- Accurate blood oxygen saturation and pulse rate.
- Anti-ambient light interference and anti- motion technology.
- Low power consumption.
- Battery display and low Power indicator Auto shutdown if no signal.

ACER BIOMEDICALS LLP

Office: Industrial Unit-80, Block-G, Sector-63, Noida, UP - 201301 (India)
Ph.: +91-120-4166921, **E-mail:** info@acerbiomed.com
Web.: www.acerbiomedicals.com, www.acerbiomed.com

ServiceHelpline
18008894232

